

# DOSECONNEC<sup>2</sup>T

by ec<sup>2</sup> Software Solutions

## TRUSTED BARCODE TECHNOLOGY

Tracking every dose from order to administration for over 4 years.

Barcode technology is a crucial tool in safely tracking medications. It enables you to ensure that the right patient is receiving the right radiopharmaceutical.



DoseConnec<sup>2</sup>t was the first in the industry and has had the ability to track radiopharmaceuticals from the order to administration for over four years.

Prior to the injection the dose can be scanned to ensure that the correct radiopharmaceutical is administered to the patient, utilizing the existing system within the hospital.

## WHY USE DOSECONNEC<sup>2</sup>T?

- **Ensure FDA Compliance with Future Mandates**
- **Reduce Administration Errors**
- **No Additional Cost**
- **Improves Prescription to Dose Accuracy**



## FUTURE DEVELOPMENTS

ec<sup>2</sup> Software Solutions is working on the next generation of barcoding technology, which will allow you to confirm both the radiopharmaceutical and dosage. This contributes further to patient safety and will continue to make ec<sup>2</sup> Software Solutions the industry leader in medical software. It's time to be proactive in radiopharmaceutical tracking.


**For more information, please contact your local Triad Isotopes pharmacy.**

# IMPLEMENTATION OF DOSECONNECT<sup>2</sup>T

**YOUR FACILITY CAN EASILY IMPLEMENT DOSECONNECT<sup>2</sup>T BY TAKING THE FOLLOWING STEPS:**



Contact your Triad Isotopes sales director to request DoseConnect<sup>2</sup>t barcoding information.



Enter the radiopharmaceutical codes provided by Triad Isotopes into your EHR/eRX database where it will be mapped to the correct radiopharmaceutical.



Scan your prescription to verify that the patient received the correct radiopharmaceutical.

**Easy as 1-2-3**

No programming required.



**Triad Isotopes**<sup>®</sup>

patient focused. | community driven.